CFX Real Time PCR system

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
The 'Plate Loading Guide' dialog box is shown with the following fields and annotations:

- Select Fluorophores...**: 選取實驗使用的螢光試劑
- Sample Type**: Unknown (dropdown) - 選取樣品種類(unknown, standard, NTC, PC, NC, NRT)
- Target Names**:
 - Load ☒ FAM <none> (dropdown) +
 - Load ☒ HEX <none> (dropdown) +
 - Load ☒ Texas <none> (dropdown) +
 - Load ☒ Cy5 <none> (dropdown) +
 - Load ☒ Quasar <none> (dropdown) +
 選取此次實驗所用的螢光試劑，並註記 target name 名稱
- Sample Names**: Load ☐ <none> (dropdown) + - 輸入樣品名稱
- Biological Group**: Load ☐ <none> (dropdown) + - 輸入組別名稱
- Replicate #**: Load ☐ 1 (dropdown) - 輸入重複性樣品設定(Horizontal 水平重複/Vertical 垂直重複)
- Technical Replicates**: Load ☐ 1 (dropdown) - 輸入重複性樣品設定(Horizontal 水平重複/Vertical 垂直重複)
- Experiment Settings...**: 設定 ref gene 與 control sample (圖十)
- Clear Replicate #**: 清除所選的重複
- Clear Wells**: 清除所選的輸入

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選擇 Reference gene

選擇 Control Sample

The image shows two side-by-side screenshots of the 'Experiment Settings' dialog box. The left window is titled 'Experiment Settings' and has a 'Targets' tab selected. It contains a table with columns: Name, Full Name, Reference, and Select To Remove. The table has four rows: 1 (18s), 2 (actin), 3 (GAPDH), and 4 (tubulin). The 'Reference' column for row 3 is checked. Below the table are buttons for 'New', 'Add', and 'Remove checked item(s)', and checkboxes for 'Show Analysis Settings' and 'Exclude the following sample types from Gene Expression analysis: NTC, NRT, Negative Control, Positive Control, Standard'. The right window is also titled 'Experiment Settings' and has a 'Targets' tab selected. It contains a table with columns: Name, Full Name, Control, and Select To Remove. The table has four rows: 1 (1), 2 (2), 3 (3), and 4 (4). The 'Control' column for row 1 is checked. Below the table are buttons for 'New', 'Add', and 'Remove checked item(s)', and checkboxes for 'Show Analysis Settings' and 'Exclude the following sample types from Gene Expression analysis: NTC, NRT, Negative Control, Positive Control, Standard'.

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The image shows a 'Starting Concentration' dialog box. It has the following fields and controls: 'Starting Concentration' (1.00E+06), 'Replicates from' (1), 'to' (3), 'Dilution Factor' (10.000), and radio buttons for 'Increasing' and 'Decreasing' (selected). Below these are 'Cancel' and 'Apply' buttons. Red arrows point from Chinese labels to the fields: '標準品的起始濃度' points to 'Starting Concentration', '序列稀釋次數' points to 'Replicates from' and 'to', and '稀釋倍數' points to 'Dilution Factor'.

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