UVP UVsolo touch Imaging System

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
The screenshot shows the software interface with the following elements:

- Gallery Tab:** A red box highlights the template selection area containing: UV-100 ms, UV-500 ms, UV-1 sec, UV-Auto Capture, White, and Manual Input.
- Capture Mode Panel:**
 - Manual Exposure:** A dropdown menu.
 - Exposure Controls:** Four vertical sliders for Min, Sec, Ms, and Micro. The Ms slider is set to 39. A red arrow points to these controls with the label "曝光時間設定".
 - Number of frames:** A control set to 1.
 - Done:** A checkmark icon at the bottom left.
- Preview Window:** A large black area showing a circular scale and a grid of horizontal lines.
- Bottom Control Bar:**
 - Stop Preview:** A button with a plus icon in a box, highlighted with a red box and a red arrow. Below it is the label "開啟/關閉預覽功能".
 - Start Capture:** A button with a camera icon in a box, highlighted with a red box and a red arrow. Below it is the label "開始拍攝".
 - Status:** "Preview Started... Exposing: 0.039 Sec".

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The screenshot displays the UVP software interface. At the top, there are 'Acquisition' and 'Gallery' tabs. Below the tabs is a gallery of three images, with the third image selected. The main window shows a large image of a clock and a ruler. On the left side, there is a list of image processing options:

Icon	Option	Chinese Label
	Histogram	灰階比
	Time Stamp	影像時間
	Saturation	過曝警示
	Pseudocolor	套色
	Invert	色差反轉

At the bottom right of the main image, there is a timestamp: Nov/17/2015 17:28:25. Below the main image, there are icons for 'Save' (存檔) and 'Print' (列印).

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

Technical Data

Specifications/features	UVP UVsolo touch
Camera/Lens	GelCam 315, 5.0 MP, with 8-48mm, f/1.2 manual zoom lens
Darkroom features	<ul style="list-style-type: none"> ▪ 11.6" large articulating touch screen ▪ Wide front door with UV safety switch ▪ Gel viewer window ▪ Epi white light ▪ Side access doors ▪ USB ports ▪ Emission filter tray <ul style="list-style-type: none"> ▪ Support for 21 CFR Part 11
Accessories included	<ul style="list-style-type: none"> ▪ Ethidium bromide filter ▪ Choice of transilluminator (Filter sizes: 20 x 20 cm, 21 x 26 cm or 25 x 26 cm; Wavelength: 1 UV 302 nm, 2UV 302/365 nm or 3 UV 254/302/365 nm, 100-115 V/60 Hz or 230 V/50 Hz)
Dimensions: H x W x D	780 x 361 x 338 mm (30.7 x 14.2 x 13.3 inch)
Wireless network capability	Wireless network capable, Wi-Fi, accessory for wired-to-ethernet connection available

Thanks For Your Attention !